

New DNA test will improve tracking of *Salmonella* food-poisoning outbreaks

28 April 2020



"*Salmonella* in a clinical or food sample may be present in very small amounts and thus requires very sensitive methods to detect them. Multiple cross-displacement amplification (MCDA) is a method that can detect very small amounts of DNA rapidly and is also performed at a single constant temperature, in contrast to the cycling of temperatures required in other methods such as PCR. This makes it a good fit for a simple, rapid, and sensitive bacterial detection test. Although an MCDA test for any *Salmonella* already exists, it does not distinguish between different serotypes," explained Professor Ruiting Lan, Ph.D., of the School of Biotechnology and Biomolecular Sciences at the University of New South Wales, Sydney, NSW, Australia.

Salmonella are gram-negative, rod-shaped, facultative anaerobic bacteria. Zhang et al developed sensitive and specific assays to detect different serotypes of *Salmonella*, paving the way for rapid serotyping directly from specimens. Credit: Public Health Image Library, US Centers for Disease Control and Prevention, James Archer (2019).

Researchers report the development of a sensitive and specific assay to detect different serotypes of *Salmonella*, paving the way for rapid serotyping directly from specimens. This improvement upon current testing methods can play a critical role in quickly tracing the origin of the infection. The report appears in *The Journal of Molecular Diagnostics*.

Recent data from the Centers for Disease Control and Prevention (CDC) indicate that food poisoning caused by *Salmonella* bacteria leads to 1.35 million infections, 26,500 hospitalizations, and 420 deaths in the United States per year. During an outbreak, the speed and simplicity of a test to detect specific types of bacteria are important for public health investigators to track down the source.

The investigators developed an MCDA assay for each of the seven serovar (subtype)-specific targets of *Salmonella*. All of these assays accurately detect as few as 10 copies of DNA and can produce results in approximately eight minutes. Importantly, these assays do not require specialized detection equipment, simplifying any future application in clinical or industrial settings. By combining these seven serovar-specific assays with the existing species assay, *Salmonella* detection can be simplified and accelerated.

"The assays developed in this study are unique because the gene markers used were selected based on analyzing thousands of genomes. Thus, these markers future proof *Salmonella* serotyping in the era of culture-independent diagnostic testing," commented Professor Lan.

Traditional methods to distinguish *Salmonella* serotypes involve growing the bacteria from samples and then testing them to assign them to a serovar. The MCDA test is faster because it does not require first growing the bacteria in culture. Rather, it can detect very small amounts of *Salmonella* DNA.

Although there are hundreds of *Salmonella* serovars, the authors chose the five most commonly occurring in Australia, which cause more than 85 percent of *Salmonella* infections in that country. However, at least two of these serovars are also the top *Salmonella* serovars throughout the world, therefore the researchers believe the results are applicable to other [geographic regions](#).

Salmonella bacteria typically live in animal and human intestines and are shed through feces. Humans may become infected by ingesting contaminated water or food such as raw or undercooked meat, poultry, eggs, or egg products. Possible signs and symptoms of *Salmonella* poisoning include nausea, vomiting, abdominal cramps, diarrhea, fever, and headache.

More information: *Journal of Molecular Diagnostics* (2020). DOI: [10.1016/j.jmoldx.2020.02.006](https://doi.org/10.1016/j.jmoldx.2020.02.006)

Provided by Elsevier

APA citation: New DNA test will improve tracking of Salmonella food-poisoning outbreaks (2020, April 28) retrieved 30 August 2022 from <https://medicalxpress.com/news/2020-04-dna-tracking-salmonella-food-poisoning-outbreaks.html>

This document is subject to copyright. Apart from any fair dealing for the purpose of private study or research, no part may be reproduced without the written permission. The content is provided for information purposes only.